



# MxDC Developer Documentation

*Release 4.1.6-gedc482b*

March 19, 2012



# CONTENTS

<b>I</b>	<b>Contents</b>	<b>1</b>
<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Prerequisites . . . . .	3
1.2	Installation . . . . .	5
1.3	Usage . . . . .	5
<b>2</b>	<b>API documentation</b>	<b>7</b>
2.1	Interfaces . . . . .	7
2.2	Devices . . . . .	8
2.3	Engines . . . . .	25
2.4	Beamline . . . . .	27
<b>3</b>	<b>Additional Information</b>	<b>29</b>
	<b>Index</b>	<b>31</b>



Part I  
CONTENTS



# INTRODUCTION

MxDC (Macromolecular Crystallography Data Collector) is a modular software package layered above the low-level beamline instrument control system. The first layer, also called the beamline control module (BCM), provides a functional abstraction of all hardware devices usually found at MX beamlines, while the graphical userface provides an integrated user-friendly interface for interactive data acquisition.

This is the documentation for MxDC and the BCM. It documents the application programming interfaces for developers wishing extend or customize the software. Some information about dependences and installation procedures is also provided.

This document is targeted at beamline administrators or software developers rather than users of MxDC. If you are a user looking for user documentation, please consult the [MxDC user guide](#)

## 1.1 Prerequisites

MxDC and the BCM have been designed and tested to work on Linux systems. Although it would probably work after considerable effort on other Unix-like system or even on Windows, we do not recommend Non-Linux systems.

MxDC and the BCM needs at least **Python 2.4** to run. It does not yet run under Python 3+. It is a pure python package which requires no compilation. Additionally, several python modules are required, without which MxDC and the BCM will not run or function correctly.

### Required modules:

- matplotlib >= 0.99
- numpy >= 1.2
- scipy >= 0.6
- python-ctypes (for python < 2.5 only)

- python-simplejson (for python < 2.5 only)
- python-imaging >= 1.1.6
- python-zope-interface >= 3.3
- Twisted >= 8.2
- avahi-tools (provides python avahi module needed by mdns module)
- dbus-python, usually installed by default on most recent distributions
- Others: (pygtk usually included with distribution)

---

**Note:** If using Twisted >= 9.0, pyasn1, and python-crypto are also required

---

### General tips:

- To avoid dependency nightmares, always prefer packages already available for your distribution. Only build where necessary.
- For installation at the CMCF beamlines, dependencies not available as part of the distribution have been built and are available in /cmcf\_apps/deps. Select rpm folder matching your architecture (x86 or x86\_64) and install the packages.
- On Scientific Linux/RHEL/CENTOS 5+, most of the packages are available from distribution repositories. Make sure you have EPEL and RPMFusion repositories enabled from the start. That way, doing yum install <package-name> should get you most of the above packages “for free”.
- A few of the packages may need to be built. It is recommended to build rpm’s and then install those, so you can easily uninstall them or update them, and it is easy to tell what versions are installed. Therefore, it is not recommended to do:

```
python setup.py install
```

Instead, do:

```
python setup.py bdist_rpm
```

which should generate an rpm package within the dist subdirectory. You can then install the rpm package normally.

### Other Programs/Libraries:

- CBFLib: this library must be installed to allow MxDC to recognize CBF diffraction image files. If it is not available, only MarCCD and SMV formatted image files will be supported.
- XREC is used for automated crystal centering and should be available in the execution path. Without it, the automatec “crystal” and “loop” centering will fail.



- CHOOCH is required for MAD and SAD scans.
- MxLIVE and Autoprocess servers are automatically discovered at run-time but can be specified in the configuration file manually if not part of the local network. Without these, automated data processing and uploading of datasets or processing reports will not work.

## 1.2 Installation

1. Uncompress the tarball archive into a target directory
2. Edit `deploy/bcm.csh` and `deploy/bcm.csh` to suit the installation
3. Create a configuration file in `etc/xxx.py` to match your beamline hardware. It is better to copy an existing one and modify it.
4. Make sure `bcm.csh` is sourced on login.

## 1.3 Usage

To run the programs:

- `mxdc` will launch the application.
- `sim-mxdc` will launch the application with simulated devices.
- `sampleviewer` will launch the hutch sample video screen only.
- `hutchviewer` will launch the hutch control screen only.
- `blconsole` will launch a beamline python console for interactive scanning and manipulation etc.



# API DOCUMENTATION

This section is automatically generated from the code.

## 2.1 Interfaces

```
class bcm.device.interfaces.IDevice(name, bases=(), attrs=None, __doc__=None,
                                     __module__=None)
```

Bases: `zope.interface.Interface`

A generic device interface.

[illegible]Bases: `bcm.device.interfaces.IDevice`

A sample automounter device interface.

[illegible]Bases: `bcm.device.interfaces.IDevice`

An X-ray beam collimator device object.

```
class bcm.device.interfaces.ICounter(name, bases=(), attrs=None, __doc__=None,  
                                     __module__=None)
```

Bases: `bcm.device.interfaces.IDevice`

An integrating counter object.

[illegible]Bases: `bcm.device.interfaces.IDevice`

A goniometer device object.

```
class bcm.device.interfaces.IImagingDetector(name, bases=(), attrs=None,
__doc__=None, __module__=None)
```

Bases: bcm.device.interfaces.IDevice

An imaging detector device for acquiring image frames.

## 2.2 Devices

### 2.2.1 Base Device

```
class bcm.device.base.BaseDevice
```

Bases: gobject.\_gobject.GObject

A generic device object class. All devices should be derived from this class. Objects of this class are instances of *gobject.GObject*.

#### Attributes:

- *pending\_devs*: a list of inactive child devices.
- *health\_manager*: A HealthManager object.
- *state\_info*: A dict containing state information.
- *name*: the name of the device

#### Signals:

- *active*: emitted when the state of the device changes from inactive to active and vice-versa. Passes a single boolean value.
- *busy*: emitted to notify listeners of a change in the busy state. A single boolean parameter is passed along with the signal. *True* means busy, *False* means not busy.
- *health*: signals an device sanity/error condition. Passes two parameters, an integer error code and a string description. The integer error codes are:
  - 0 : No error
  - 1 : MINOR, no impact to device functionality.
  - 2 : MARGINAL, no immediate impact. Attention may soon be needed.
  - 4 : SERIOUS, functionality impacted but recovery is possible.
  - 8 : CRITICAL, functionality broken, recovery is not possible.
  - 16 : DISABLED, device has been manually disabled.

- *message*: signal for sending messages from the device to any listeners. Messages are passed as a single string parameter.

### Dynamic Attributes:

- *<signal>\_state*: for each signal defined in derived classes, the last transmitted data corresponding to the current state can be obtained through the dynamic variable with the suffix “\_state” added to the signal name. For example:

```
>>> obj.active_state
>>> True
```

### Dynamic Methods:

- *is\_<signal>()*: for each boolean signal defined in derived classes, the current state can be obtained by adding the prefix “**is\_**” to the signal name. For example:

```
>>> obj.is_active()
>>> True
```

`get_state()`

Obtain a copy of the device state.

**Returns:** A dict mapping state keys to their corresponding values. Entries contain at least the following entries:

- *active* : Boolean
- *busy*: Boolean
- *health*: tuple(int, string)
- *message*: string

`set_state(**kwargs)`

Set the state of the device and emit the corresponding signal.

**Kwargs:** Keyworded arguments follow the same conventions as the state dictionary and correspond to any signals defined for the device. Signals must be previously defined as supported signals of the device.

For example:

```
mydevice.set_state(active=True, busy=False,
                   health=(1, 'error', 'too hot'),
                   message="the device is overheating")
```

`add_pv(*args, **kwargs)`

Add a process variable (PV) to the device.

Create a new process variable, add it to the device. Keyworded arguments should be the same as those expected for instantiating a `bcm.protocol.ca.PV` object.

**Returns:** A reference to the created `bcm.protocol.ca.PV` object.

`add_devices(*devices)`

Add one or more devices as children of this device.

**class** `bcm.device.base.HealthManager(**kwargs)`

Bases: object

Manages the health states. The object enables registration and removal of error states and consistent reporting of health based on all currently active health issues.

Takes key worded arguments. The keyword name is the context, and the value is an error string to be returned instead of the context name with all health information for the given context.

`register_messages(**kwargs)`

Update or add entries to the context message register

`add(severity, context, msg=None)`

Adds an error state of the given context as a string and a severity value as a integer. If a message is given, it will be stored and used instead of the context name. Only one message per context type is allowed. Use a different context if you want different messages.

`remove(context)`

Remove all errors of the given context as a string.

`get_health()`

Generate an error code and string based on all the currently registered errors within the health registry.

## 2.2.2 Automounter

**class** `bcm.device.automounter.BasicAutomounter`

Bases: `bcm.device.base.BaseDevice`

Basic Automounter object.

An automounter object contains a number of AutomounterContainers.

**Signals:**

- *state*: transmits changes in automounter states. The signal data is a dictionary with three entries {'busy': <boolean>,
- *enabled*: <boolean>, 'needs': [<str>, <str>, ...]}.

- *mounted*: emitted when a sample is mounted or dismounted. The data transmitted is a tuple of the form (<port no.>, <barcode>) when mounting and *None* when dismounting.
- **progress: notifies listeners of automounter progress. Data is a tuple of the form (<% complete>, <description>).**

`abort()`

Abort all actions.

`probe()`

Command the automounter the probe for containers and port states.

`mount(port, wash=False, wait=False)`

Mount the sample at the specified port, optionally washing the sample in the process. Does nothing if the requested port is already mounted. Dismounts the mounted sample prior to mounting if a another sample is already mounted.

**Args:**

- *port* (str): Address to mount.

**Kwargs:**

- *wash* (bool): Whether to wash or not (default is False)
- *wait* (bool): Run asynchronously or block (default is async, False)

**Returns:** bool. True if the requested sample is successfully mounted, and False otherwise.

`dismount(port=None, wait=False)`

Dismount the sample into the specified port if provided, otherwise dismount it to the original port from which it was mounted.

**Kwargs:**

- *port* (str): Destination address to dismount to, default is original port
- *wait* (bool): Run asynchronously or block (default is async, False)

`wait(start=True, stop=True, timeout=240)`

Wait for the automounter

**Kwargs:**

- *start* (bool): Wait for automounter to start.
- *stop* (bool): Wait for automounter to stop.
- *timeout* (int): Maximum time to wait

**Returns:** True if the wait was successful. False if the wait timed-out.

`is_mountable(port)`

Check if a sample location can be mounted safely. Does not guarantee that the port is actually mountable, only that the the port was marked as mountable during the last probe operation.

**Args:** *port* (str): The sample location to check.

**Returns:** True or False.

`is_mounted(port=None)`

Check if any sample or a specific sample location is currently mounted.

**Kwargs:** *port* (str): The sample location to check.

**Returns:** True if a port is specified and is mounted or if no port is specified and any sample is mounted. False otherwise.

`get_port_state(port)`

Obtain the detailed state of the specified sample location.

**Args:** *port* (str): The sample location

**Returns:** int. one of:

- 0 -- PORT\_EMPTY
- 1 -- PORT\_GOOD
- 2 -- PORT\_UNKNOWN
- 3 -- PORT\_MOUNTED
- 4 -- PORT\_JAMMED
- 5 -- PORT\_NONE

### 2.2.3 Counters

**class** `bcm.device.counter.Counter(pv_name, zero=0.0)`

Bases: `bcm.device.base.BaseDevice`

EPICS based Counter Device objects. Enables counting and averaging of process variables over given time periods.

**Args:** *pv\_name* (str): Process Variable name.

**Kwargs:** *zero* (float): Zero offset value. Defaults to 0.0;

**Returns** float. The average process variable value during the count time.

`count(t)`

Count for the specified amount of time and return the numeric value corresponding to the average count. This method blocks for the specified time.

**Args:** *t* (float): averaging time in seconds.



**Returns** float. The average process variable value during the count time.

## 2.2.4 Detector

```
class bcm.device.detector.MXCCDImager(name, size, resolution, detector_type='MX300')
```

Bases: bcm.device.base.BaseDevice

MX Detector object for EPICS based Rayonix CCD detectors at the CLS.

**Args:** *name* (str): Root name of the EPICS record Process Variables. *size* (int): The size in pixels of the detector. Assumes square detectors. *resolution* (float): The pixel size in

**Kwargs:** *detector\_type* (str): The type of detector. e.g. "MX300" for Rayonix MX CCD 300.

```
initialize(wait=True)
```

Initialize the detector and take background images if necessary. This method does not do anything if the device is already initialized.

**Kwargs:** *wait* (bool): If true, the call will block until initialization is complete.

```
start(first=False)
```

Start acquiring.

**Kwargs:** *first* (bool): Specifies whether this is the first of a series of acquisitions. This is used to customize the behaviour for the first.

```
stop()
```

Stop and Abort the current acquisition.

```
save(wait=False)
```

Save the current buffers according to the current parameters.

**Kwargs:** *wait* (bool): If true, the call will block until the save operation is complete.

```
get_origin()
```

Obtain the detector origin/beam position in pixels.

**Returns:** tuple(x, y) corresponding to the beam-x and beam-y coordinates.

```
wait(state='idle')
```

Wait until the detector reaches a given state.

**Kwargs:** *state* (str): The state to wait for. Default 'idle'.

```
set_parameters(data)
```

Set the detector parameters for the image header and file names.

**Args:** *data* (dict): A dictionary of key value pairs for the parameters. supported parameters are:

- *filename* (str), Output file name of the image.
- *directory* (str), Directory name to store image.
- *beam\_x* (int), Detector X-origin in pixels.
- *beam\_y* (int), Detector Y-origin in pixels.
- *distance* (float), Detector distance in mm.
- *exposure\_time* , Exposure time in seconds.
- *axis* (str), Spindle rotation axis.
- *wavelength* (float), Wavelength of radiation in Angstroms.
- *delta\_angle* (float), Delta oscillation angle in deg.
- *frame\_number* (int), Frame number.
- *name* (str), File name prefix for the image.
- *start\_angle* (float), Starting spindle position of image in deg.
- *energy* (float), Wavelength of radiation in KeV.
- *comments* (str), File comments.

### 2.2.5 Diffractometer

**class** bcm.device.diffractometer.Diffractometer(*distance*, *two\_theta*,  
name='Diffractometer')

Bases: bcm.device.base.BaseDevice

A Container device object which groups a detector distance and detector swing-out.

**Args:**

- *distance* (class::interfaces.IMotor provider): device which controls the detector distance.
- *two\_theta* (class::interfaces.IMotor provider): device which controls the detector swing-out angle.

**Kwargs:** *name* (str): The name of the device group.

wait()

Wait for both child devices to finish moving.

stop()

Stop both child devices.

## 2.2.6 Goniometer

**class** bcm.device.goniometer.BackLight(*name*)

Bases: bcm.device.misc.BasicShutter

A specialized in-out actuator for pneumatic OAV backlight at the CLS.

**Args:** *-name* (str): Root PV name of EPICS record.

**class** bcm.device.goniometer.Goniometer(*name, blname, mnt\_cmd, minibeam*)

Bases: bcm.device.goniometer.GoniometerBase

EPICS based Parker-type Goniometer at the CLS 08ID-1.

**Args:**

- *name* (str): PV name of goniometer EPICS record.
- *blname* (str): PV name for Backlight PV.
- *mnt\_cmd* (str): PV name for toggling mount mode.
- *minibeam* (str): PV name for minibeam motor.

configure(*\*\*kwargs*)

Configure the goniometer to perform an oscillation scan.

**Kwargs:** *time* (float): Exposure time in seconds *delta* (float): Delta oscillation range in deg *angle* (float): Starting angle of oscillation in deg

set\_mode(*mode, wait=False*)

Set the mode of the goniometer environment.

**Args:**

- *mode* (str) one of:
  - “CENTERING” : Prepare for centering
  - “MOUNTING” : Prepare for mounting/dismounting samples
  - “COLLECT” : Prepare for data collection
  - “BEAM” : Inspect the beam
  - “SCANNING” : Prepare for scanning and fluorescence measurements

**Kwargs:**

- *wait* (bool): if True, block until the mode is completely changed.

scan(*wait=True*)

Perform an oscillation scan according to the currently set parameters

**Kwargs:**

- *wait* (bool): if True, wait until the scan is complete otherwise run asynchronously.

**class** bcm.device.goniometer.MD2Goniometer(*name*)

Bases: bcm.device.goniometer.GoniometerBase

EPICS based MD2-type Goniometer at the CLS 08B1-1.

**Args:**

- *name* (str): Root PV name of the goniometer EPICS record.

configure(\*\**kwargs*)

Configure the goniometer to perform an oscillation scan.

**Kwargs:**

- *time* (float): Exposure time in seconds
- *delta* (float): Delta oscillation range in deg
- *angle* (float): Starting angle of oscillation in deg

set\_mode(*mode*, *wait=False*)

Set the mode of the goniometer environment.

**Args:**

- *mode* (str) one of:
  - “CENTERING” : Prepare for centering
  - “MOUNTING” : Prepare for mounting/dismounting samples
  - “COLLECT” : Prepare for data collection
  - “BEAM” : Inspect the beam
  - “SCANNING” : Prepare for scanning and fluorescence measurements

**Kwargs:**

- *wait* (bool): if True, block until the mode is completely changed.

scan(*wait=True*)

Perform an oscillation scan according to the currently set parameters

**Kwargs:**

- *wait* (bool): if True, wait until the scan is complete otherwise run asynchronously.

stop()

Stop and abort the current scan if any.

## 2.2.7 Multi-Channel Analyzer

**class** bcm.device.mca.XFlashMCA(*name*, *nozzle*=None, *channels*=4096)

Bases: bcm.device.mca.BasicMCA

mcaRecord based single element fluorescence detector object.

**Args:**

- *name* (str): Root PV name of the mcaRecord.

**Kwargs:**

- *nozzle* (bcm.device.misc.Positioner): Nozzle positioner.
- *channels* (int): Number of channels.

configure(\*\**kwargs*)

Configure the detector for data acquisition.

**Kwargs:**

- *retract* (bool): True means retract the nozzle.
- *roi* (tuple(int, int)): bounding channels enclosing region of interest.
- *energy* (float): an energy value in keV around which to construct a region of interest. The ROI is calculated as a 150 eV window around the requested energy. If both *roi* and *energy* are given, *energy* takes precedence.
- *cooling* (bool): True means enable detector cooling, False means disable it.

**class** bcm.device.mca.VortexMCA(*name*, *channels*=2048)

Bases: bcm.device.mca.BasicMCA

EPICS based 4-element Vortex ME4 detector object.

**Args:** *name* (str): Root PV name of EPICS record.

**Kwargs:** *channels* (int): Number of channels.

**class** bcm.device.mca.SimMultiChannelAnalyzer(*name*, *channels*=4096)

Bases: bcm.device.mca.BasicMCA

Simulated single channel MCA detector.

**Args:** *name* (str): Name of device.

**Kwargs:** *channels* (int): Number of channels.

### 2.2.8 Monochromator

```
class bcm.device.monochromator.Monochromator(simple_energy, energy,
                                             mostab=None)
    Bases: bcm.device.base.BaseDevice
    wait()
    stop()
```

### 2.2.9 Motors

**exception** `bcm.device.motor.MotorError`  
 Bases: `exceptions.Exception`  
 Base class for errors in the motor module.

```
class bcm.device.motor.MotorBase(name)
    Bases: bcm.device.base.BaseDevice

    Base class for motors.
```

**Signals:**

- *changed* (float): Emitted everytime the position of the motor changes. Data contains the current position of the motor.
- *timed-change* (tuple(float, float)): Emitted everytime the motor changes. Data is a 2-tuple with the current position and the timestamp of the last change.

```
class bcm.device.motor.Motor(pv_name, motor_type)
    Bases: bcm.device.motor.MotorBase
```

Motor object for EPICS based motor records.

**Args:**

- *pv\_name* (str): Root PV name for the EPICS record.
- *motor\_type* (str): Type of EPICS motor record. Accepted values are:
  - "vme" - CLS VME58 and MaxV motor record without encoder support.
  - "vmeenc" - CLS VME58 and MaxV motor record with encoder support.
  - "cls" - OLD CLS motor record.
  - "pseudo" - CLS PseudoMotor record.

`get_position()`  
Obtain the current position of the motor in device units.

**Returns:** float.

`configure(**kwargs)`

Configure the motor. Currently only allows changing the calibration flag and resetting the position for “vme” and “vmeenc” motors only.

**Kwargs:**

- *calib* (int): 0 removes the calibration flag, anything else sets it.
- *reset* (float): position to recalibrate the current motor to.

`move_to(pos, wait=False, force=False)`

Request the motor to move to an absolute position. By default, the command will not be sent to the motor if its current position is the same as the requested position within its preset precision. In addition the command will not be sent if the motor health severity is not zero (GOOD).

**Args:**

- *pos* (float): Target position to move to.

**Kwargs:**

- *wait* (bool): Whether to wait for move to complete or return immediately.
- *force* (bool): Force a command to be sent to the motor even if the target is the same as the current position.

`move_by(val, wait=False, force=False)`

Similar to `move_to()`, except request the motor to move by a relative amount.

**Args:**

- *val* (float): amount to move position by.

**Kwargs:**

- *wait* (bool): Whether to wait for move to complete or return immediately.
- *force* (bool): Force a command to be sent to the motor even if the target is the same as the current position.

`stop()`

Stop the motor from moving.

`wait(start=True, stop=True)`

Wait for the motor busy state to change.

**Kwargs:**

- *start* (bool): Wait for the motor to start moving.
- *stop* (bool): Wait for the motor to stop moving.

`get_settings()`

Obtain the motor settings for saving/restore purposes.

**Returns:** A dict.

**class** `bcm.device.motor.VMEMotor(name)`

Bases: `bcm.device.motor.Motor`

Convenience class for “vme” type motors.

**Args:**

- *name* (str): Root PV name of the motor record.

**class** `bcm.device.motor.ENCMotor(name)`

Bases: `bcm.device.motor.Motor`

Convenience class for “vmeenc” type motors.

**Args:**

- *name* (str): Root PV name of the motor record.

**class** `bcm.device.motor.CLSMotor(name)`

Bases: `bcm.device.motor.Motor`

Convenience class for “cls” type motors.

**Args:**

- *name* (str): Root PV name of the motor record.

**class** `bcm.device.motor.PseudoMotor(name)`

Bases: `bcm.device.motor.Motor`

Convenience class for “pseudo” type motors.

**Args:**

- *name* (str): Root PV name of the motor record.

**class** `bcm.device.motor.BraggEnergyMotor(name, enc=None, motor_type='vme')`

Bases: `bcm.device.motor.Motor`

A specialized energy motor for using just the monochromator bragg angle.

**Args:**

- *name* (str): Root PV name of motor record.

**Kwargs:**

- *enc* (str): PV name for an optional encoder feedback value from which to read the energy value.
- *motor\_type* (str): Type of EPICS motor record. Accepted values are:



"vme" - CLS VME58 and MaxV motor record without encoder support.  
"vmeenc" - CLS VME58 and MaxV motor record with encoder support.  
"cls" - OLD CLS motor record.  
"pseudo" - CLS PseudoMotor record.

**class** bcm.device.motor.FixedLine2Motor(*x, y, slope, intercept, linked=False*)

Bases: bcm.device.motor.MotorBase

A specialized fixed offset pseudo-motor for moving two motors along a straight line.

**Args:**

- *x* (MotorBase): x-axis motor.
- *y* (MotorBase): y-axis motor.
- *slope* (float): slope of the line.
- *intercept* (float): y-intercept of the line.

**Kwargs:**

- *linked* (bool): Whether the two motors are linked. Two motors are linked if they can not be moved at the same time.

get\_position()

Obtain the position of the *x* motor only.

**class** bcm.device.motor.RelVerticalMotor(*y1, y2, omega, offset=0.0*)

Bases: bcm.device.motor.MotorBase

A specialized pseudo-motor for moving an x-y stage attached to a rotating axis vertically. Such as a centering table attached to a goniometer. The current position is always zero and all moves are relative.

**Args:**

- *y1* (MotorBase): The first motor which is moves vertically when the angle of the axis is at zero.
- *y2* (MotorBase): The second motor which is moves horizontally when the angle of the axis is at zero.
- *omega* (MotorBase): The motor for the rotation axis.

**Kwargs:**

- *offset* (float): An angle correction to apply to the rotation axis position to make *y1* vertical.

## 2.2.10 Video

**class** bcm.device.video.VideoSrc(*name*='Basic Camera', *maxfps*=5.0)

Bases: bcm.device.base.BaseDevice

Base class for all Video Sources. Maintains a list of listeners (sinks) and updates each one when the video frame changes.

Kwargs: *name* (str): Camera name. *maxfps* (float): Maximum frame rate.

*add\_sink(sink)*

Add a video sink.

**Args:** *sink* (bcm.device.interfaces.IVideoSink provider).

*del\_sink(sink)*

Remove a video sink.

**Args:** *sink* (bcm.device.interfaces.IVideoSink provider).

*start()*

Start producing video frames.

*stop()*

Start producing video frames.

*get\_frame()*

Obtain the most recent video frame.

**Returns:** A Image.Image (Python Imaging Library) image object.

## 2.2.11 Miscellaneous

**class** bcm.device.misc.PositionerBase

Bases: bcm.device.base.BaseDevice

Base class for a simple positioning device.

**Signals:**

- *changed* : Data is the new value of the device.

*set(value)*

Set the value. Args:

- *value* : New value to set.

*get()*

**Returns:**

- The current value.

**class** bcm.device.misc.Positioner(*name*, *fbk\_name*=None, *scale*=100, *units*='')

Bases: bcm.device.misc.PositionerBase

Simple EPICS based positioning device.

**Args:**

- *name* (str): Name of positioner PV for setting the value

**Kwargs:**

- *fbk\_name* (str): Name of PV for getting current value. If not provided, the same PV will be used to both set and get.
- *scale* (float): A percentage to scale the set and get values by.
- *units* (str): The units of the value.

**class** bcm.device.misc.PositionerMotor(*positioner*)

Bases: bcm.device.motor.MotorBase

Adapts a positioner so that it behaves like a Motor (ie, provides the *IMotor* interface.

**Args:**

- *positioner* (PositionerBase)

**class** bcm.device.misc.MotorShutter(*motor*)

Bases: bcm.device.base.BaseDevice

Used for CMCF1 cryojet Motor

*is\_open()*

Convenience function for open state

**class** bcm.device.misc.DiskSpaceMonitor(*descr*, *path*, *warn*=0.8, *critical*=0.9,  
*freq*=5.0)

Bases: bcm.device.base.BaseDevice

An object which periodically monitors a given path for available space.

**Args:**

- *descr* (str): A description.
- *path* (str): Path to monitor.

**Kwargs:**

- *warn* (float): Warn if fraction of used space goes above this value.
- *critical* (float): Raise an error if fraction of used space goes above this value.

*-freq* (float): Frequency in minutes to check disk usage. Default (5)

**class** bcm.device.cryojet.CryojetNozzle(*name*)  
Bases: bcm.device.misc.BasicShutter

A specialized in-out actuator for pneumatic cryoject nozzles at the CLS.

**Args:** *name* (str): Process variable root name.

**class** bcm.device.cryojet.Cryojet(*cname, lname, nname=""*)  
Bases: bcm.device.base.BaseDevice

EPICS Based cryoject device object at the CLS.

**Args:**

- *cname* (str): Root name for EPICS cryojet record.
- *lname* (str): root name for EPICS cryo-level controller record.

**Kwargs:**

- *nname* (str): Root name of the EPICS cryoject nozzle record.

stop\_flow()

Stop the flow of the cold nitrogen stream. The current setting for flow rate is saved.

resume\_flow()

Restores the flow rate to the previously saved setting.

**class** bcm.device.diagnostics.DiagnosticBase(*descr*)  
Bases: gobject.\_gobject.GObject

Base class for diagnostics.

**Signals:**

- *status* (tuple(int,str)): Changes in status. The first entry is a represents the severity while the second is a status message. Severity values can be one of:

- 0 - DIAG\_STATUS\_GOOD
- 1 - DIAG\_STATUS\_WARN
- 2 - DIAG\_STATUS\_BAD
- 3 - DIAG\_STATUS\_UNKNOWN
- 4 - DIAG\_STATUS\_DISABLED

get\_status()

Check the state of the diagnostic.

**Returns:** tuple(int, str). The string is a status description while the value of the integer corresponds to the status severity

**class** bcm.device.diagnostics.DeviceDiag(*device*, *descr=None*)

Bases: bcm.device.diagnostics.DiagnosticBase

A diagnostic object for generic devices which emits a warning when the device health is not good and an error when it is disconnected or disabled.

**Args:** *device* (a class::*device.base.BaseDevice* object) the device to monitor.

**Kwargs:** *descr* (str): Short description of the diagnostic.

**class** bcm.device.diagnostics.ServiceDiag(*service*, *descr=None*)

Bases: bcm.device.diagnostics.DiagnosticBase

A diagnostic object for generic services which emits an error when it is disconnected or disabled.

**Args:** *service* (a class::*service.base.BaseService* object) the service to monitor.

**Kwargs:** *descr* (str): Short description of the diagnostic.

## 2.3 Engines

### 2.3.1 Diffraction Engine

### 2.3.2 Spectroscopy Engine

This module defines classes and interfaces for X-Ray fluorescence.

### 2.3.3 Centering Engine

Created on Jan 26, 2010

@author: michel

### 2.3.4 Scanning Engine

**exception** bcm.engine.scanning.ScanError

Bases: exceptions.Exception

Scan Error.

**class** bcm.engine.scanning.AbsScan(*mtr*, *start\_pos*, *end\_pos*, *steps*, *cntr*, *t*, *i0=None*)

Bases: bcm.engine.scanning.BasicScan

An absolute scan of a single motor.

```
class bcm.engine.scanning.AbsScan2(mtr1, start_pos1, end_pos1, mtr2, start_pos2,  
                                   end_pos2, steps, cntr, t, i0=None)
```

Bases: bcm.engine.scanning.BasicScan

An Absolute scan of two motors.

```
class bcm.engine.scanning.RelScan(mtr, start_offset, end_offset, steps, cntr, t,  
                                  i0=None)
```

Bases: bcm.engine.scanning.AbsScan

A relative scan of a single motor.

```
class bcm.engine.scanning.CntScan(mtr, start_pos, end_pos, cntr, t, i0=None)
```

Bases: bcm.engine.scanning.BasicScan

A Continuous scan of a single motor.

```
class bcm.engine.scanning.RelScan2(mtr1, start_offset1, end_offset1, mtr2,  
                                   start_offset2, end_offset2, steps, cntr, t,  
                                   i0=None)
```

Bases: bcm.engine.scanning.AbsScan2

A relative scan of a two motors.

```
class bcm.engine.scanning.GridScan(mtr1, start_pos1, end_pos1, mtr2, start_pos2,  
                                   end_pos2, steps, cntr, t, i0=None)
```

Bases: bcm.engine.scanning.BasicScan

A absolute scan of two motors in a grid.

## 2.3.5 Fitting Engine

```
bcm.engine.fitting.peak_fit(x, y, target='gaussian')
```

**Returns the coefficients for the target function** 0 = H - Height of peak 1 = L - FWHM of Voigt profile (Wertheim et al) 2 = P - Position of peak 3 = n - lorentzian fraction, gaussian offset

Success (boolean)

```
bcm.engine.fitting.histogram_fit(x, y)
```

calcfwhm(x,y) - with input x,y vector this function calculates fwhm and returns (fwhm,xpeak,ymax, fwhm\_x\_left, fwhm\_x\_right) x - input independent variable y - input dependent variable return ymax, fwhm, xpeak, x\_hpeak[0], x\_hpeak[1], cema cema is center of mass

success boolean

## 2.3.6 Scripting Engine

**exception** `bcm.engine.scripting.ScriptError`

Bases: `exceptions.Exception`

Exceptions for Scripting Engine.

## 2.3.7 Miscellaneous Engines

Created on Oct 25, 2010

@author: michel This module defines classes for Optimizers.

**class** `bcm.engine.optimizer.PitchOptimizer(name, pitch_func, min_count=0.0)`

Bases: `bcm.device.base.BaseDevice`

Pitch Optimizer for 08B1-1

Created on Jan 26, 2010

@author: michel

## 2.4 Beamline

### 2.4.1 Interfaces

**class** `bcm.beamline.interfaces.IBeamline(name, bases=(), attrs=None,`

`__doc__=None, __module__=None)`  
Bases: `zope.interface.Interface`

A beamline object.

### 2.4.2 MX Beamline

**class** `bcm.beamline.mx.MXBeamline(console=False)`

Bases: `object`

MX Beamline(Macromolecular Crystallography Beamline) objects

Initializes a MXBeamline object from a python configuration file. The configuration file is loaded as a python module and follows the following conventions:

- Must be named the same as BCM\_BEAMLINe environment variable followed by .py and placed in the directory defined by BCM\_CONFIG\_PATH. For example if the BCM\_BEAMLINe is '08B1', the module should be '08B1.py'

- Optionally will also load a local module defined in the file \$(BCM\_BEAMLINE)\_local.py e.g '08B1\_local.py for the above example.

- **Global Variables:** BEAMLINE\_NAME = Any string preferably without spaces  
 BEAMLINE\_TYPE = Only the string 'MX' for now  
 BEAMLINE\_ENERGY\_RANGE = A tuple of 2 floats for low and hi energy limits  
 BEAMLINE\_GONIO\_POSITION = Goniometer orientation according to XREC (i.e 1,2,3 etc)  
 DEFAULT\_EXPOSURE = A float for the default exposure time  
 DEFAULT\_ATTENUATION = A float for attenuation in %  
 DEFAULT\_BEAMSTOP = Default beam-stop position  
 SAFE\_BEAMSTOP = Safe Beam-stop position during mounting  
 XRF\_BEAMSTOP = Beam-stop position for XRF scans  
 LIMS\_API\_KEY = A string  
 MISC\_SETTINGS = A dictionary containing any other key value pairs

will be available as beamline.config['misc']

**DEVICES = A dictionary mapping device names to device objects.** See SIM.py for a standard set of names.

**CONSOLE\_DEVICES = Same as above but only available in the console** in addition to the above

SERVICES = A dictionary mapping service names to service client objects  
 BEAMLINE\_SHUTTERS = A sequence of shutter device names for all shutters

required to allow beam to the end-station in the order in which they have to be opened.

Kwargs: console (bool): Whether the beamline is being used within a console or not. Used internally to register CONSOLE\_DEVICES if True. Default is False.

MXBeamline.setup()

Setup and register the beamline devices from configuration files.

### 2.4.3 Remote Client

**class** bcm.beamline.remote.BeamlineClient

Bases: gobject.\_gobject.GObject

BeamlineClient.setup(\*args, \*\*kwargs)



## ADDITIONAL INFORMATION

If you can't find the information you're looking for, have a look at the index or try to find it using the search function:

- *genindex*
- *search*



# INDEX

## A

abort() (bcm.device.automounter.BasicAutomounter  
method), 11  
AbsScan (class in bcm.engine.scanning), 25  
AbsScan2 (class in bcm.engine.scanning),  
25  
add() (bcm.device.base.HealthManager  
method), 10  
add\_devices() (bcm.device.base.BaseDevice  
method), 10  
add\_pv() (bcm.device.base.BaseDevice  
method), 9  
add\_sink() (bcm.device.video.VideoSrc  
method), 22

## B

BackLight (class in  
bcm.device.goniometer), 15  
BaseDevice (class in bcm.device.base), 8  
BasicAutomounter (class in  
bcm.device.automounter), 10  
bcm.device.automounter (module), 10  
bcm.device.base (module), 8  
bcm.device.counter (module), 12  
bcm.device.cryojet (module), 23  
bcm.device.detector (module), 13  
bcm.device.diagnostics (module), 24  
bcm.device.diffraction (module), 14  
bcm.device.goniometer (module), 15  
bcm.device.mca (module), 17  
bcm.device.misc (module), 22  
bcm.device.monochromator (module), 18  
bcm.device.motor (module), 18  
bcm.device.video (module), 22

bcm.engine.auto (module), 27  
bcm.engine.autochooch (module), 27  
bcm.engine.centering (module), 25  
bcm.engine.diffraction (module), 25  
bcm.engine.fitting (module), 26  
bcm.engine.optimizer (module), 27  
bcm.engine.scanning (module), 25  
bcm.engine.scripting (module), 27  
bcm.engine.snapshot (module), 27  
bcm.engine.spectroscopy (module), 25  
BraggEnergyMotor (class in  
bcm.device.motor), 20

## C

CLSMotor (class in bcm.device.motor), 20  
CntScan (class in bcm.engine.scanning), 26  
configure() (bcm.device.goniometer.Goniometer  
method), 15  
configure() (bcm.device.goniometer.MD2Goniometer  
method), 16  
configure() (bcm.device.mca.XFlashMCA  
method), 17  
configure() (bcm.device.motor.Motor  
method), 18  
count() (bcm.device.counter.Counter  
method), 12  
Counter (class in bcm.device.counter), 12  
Cryojet (class in bcm.device.cryojet), 24  
CryojetNozzle (class in bcm.device.cryojet),  
23

## D

del\_sink() (bcm.device.video.VideoSrc  
method), 22

DeviceDiag	(class in bcm.device.diagnostics), 24	in H
DiagnosticBase	(class in bcm.device.diagnostics), 24	HealthManager (class in bcm.device.base), 10
Diffractionmeter	(class in bcm.device.diffractionmeter), 14	histogram_fit() (in module bcm.engine.fitting), 26
DiskSpaceMonitor	(class in bcm.device.misc), 23	I
dismount() (bcm.device.automounter.BasicAutomounter method), 11		IAutomounter (class in bcm.device.interfaces), 7
E		ICollimator (class in bcm.device.interfaces), 7
ENCMotor (class in bcm.device.motor), 20		ICounter (class in bcm.device.interfaces), 7
F		IDevice (class in bcm.device.interfaces), 7
FixedLine2Motor (class in bcm.device.motor), 21	in	IGoniometer (class in bcm.device.interfaces), 7
G		IImagingDetector (class in bcm.device.interfaces), 7
get() (bcm.device.misc.PositionerBase method), 22		initialize() (bcm.device.detector.MXCCDIImager method), 13
get_frame() (bcm.device.video.VideoSrc method), 22		is_mountable() (bcm.device.automounter.BasicAutomounter method), 11
get_health() (bcm.device.base.HealthManager method), 10		is_mounted() (bcm.device.automounter.BasicAutomounter method), 12
get_origin() (bcm.device.detector.MXCCDIImager method), 13		is_open() (bcm.device.misc.MotorShutter method), 23
get_port_state() (bcm.device.automounter.BasicAutomounter method), 12		M
get_position() (bcm.device.motor.FixedLine2Motor method), 21		MD2Goniometer (class in bcm.device.goniometer), 16
get_position() (bcm.device.motor.Motor method), 18		Monochromator (class in bcm.device.monochromator), 18
get_settings() (bcm.device.motor.Motor method), 19		Motor (class in bcm.device.motor), 18
get_state() (bcm.device.base.BaseDevice method), 9		MotorBase (class in bcm.device.motor), 18
get_status() (bcm.device.diagnostics.DiagnosticBase method), 24		MotorError, 18
Goniometer (class in bcm.device.goniometer), 15	in	MotorShutter (class in bcm.device.misc), 23
GridScan (class in bcm.engine.scanning), 26		mount() (bcm.device.automounter.BasicAutomounter method), 11
		move_by() (bcm.device.motor.Motor method), 19
		move_to() (bcm.device.motor.Motor method), 19
		MXCCDIImager (class in bcm.device.detector), 13

## P

peak\_fit() (in module bcm.engine.fitting), 26

PitchOptimizer (class in bcm.engine.optimizer), 27

Positioner (class in bcm.device.misc), 22

PositionerBase (class in bcm.device.misc), 22

PositionerMotor (class in bcm.device.misc), 23

probe() (bcm.device.automounter.BasicAutomounter method), 11

PseudoMotor (class in bcm.device.motor), 20

## R

register\_messages() (bcm.device.base.HealthManager method), 10

RelScan (class in bcm.engine.scanning), 26

RelScan2 (class in bcm.engine.scanning), 26

RelVerticalMotor (class in bcm.device.motor), 21

remove() (bcm.device.base.HealthManager method), 10

resume\_flow() (bcm.device.cryojet.Cryojet method), 24

## S

save() (bcm.device.detector.MXCCDIImager method), 13

scan() (bcm.device.goniometer.Goniometer method), 15

scan() (bcm.device.goniometer.MD2Goniometer method), 16

ScanError, 25

ScriptError, 27

ServiceDiag (class in bcm.device.diagnostics), 25

set() (bcm.device.misc.PositionerBase method), 22

set\_mode() (bcm.device.goniometer.Goniometer method), 15

set\_mode() (bcm.device.goniometer.MD2Goniometer method), 16

set\_parameters() (bcm.device.detector.MXCCDIImager method), 13

set\_state() (bcm.device.base.BaseDevice method), 9

SimMultiChannelAnalyzer (class in bcm.device.mca), 17

start() (bcm.device.detector.MXCCDIImager method), 13

start() (bcm.device.video.VideoSrc method), 22

stop() (bcm.device.detector.MXCCDIImager method), 13

stop() (bcm.device.diffractionmeter.Diffractionmeter method), 14

stop() (bcm.device.goniometer.MD2Goniometer method), 16

stop() (bcm.device.monochromator.Monochromator method), 18

stop() (bcm.device.motor.Motor method), 19

stop() (bcm.device.video.VideoSrc method), 22

stop\_flow() (bcm.device.cryojet.Cryojet method), 24

## V

VideoSrc (class in bcm.device.video), 22

VMEMotor (class in bcm.device.motor), 20

VortexMCA (class in bcm.device.mca), 17

## W

wait() (bcm.device.automounter.BasicAutomounter method), 11

wait() (bcm.device.detector.MXCCDIImager method), 13

wait() (bcm.device.diffractionmeter.Diffractionmeter method), 14

wait() (bcm.device.monochromator.Monochromator method), 18

wait() (bcm.device.motor.Motor method), 19

## X

XFlashMCA (class in bcm.device.mca), 17